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# EN: SUPERNOVA by GLOV – User Manual

## Contents

1. General Device Information
2. Product Overview
3. Description of the Device
4. Before You Start
5. How to Use SUPERNOVA
6. Warnings
7. Potential Benefits
  - 7.1 Tips for Use
8. Using & Caring for SUPERNOVA
  - 8.1 General Use
  - 8.2 Cleaning
  - 8.3 Storing
  - 8.4 Troubleshooting
9. Specifications
10. Travel and International Use
11. Limited Warranty – 2 Years
12. Disposal
13. Symbols
14. Contact Information

## 1. General Device Information

SUPERNOVA LED Patch – a star-shaped beauty accessory designed for precise, at-home spot care. Combining blue light (415 nm) and red light (630 nm), it brings together two proven wavelengths in a compact and reusable format. With its fun design and portable size, Supernova makes daily light care simple and effective. Each full charge allows 8-12 quick 3-minute sessions, perfect for incorporating into your skincare routine whenever and wherever you need it.

**⚠ Note:** This is a beauty device for personal home use only. It is not a medical device and is not intended to diagnose, treat, or cure any medical condition.

## 2. Product Overview

| SUPERNOVA Case                                                                                                                                | SUPERNOVA USB                                                                                                            |
|-----------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------|
| 1 × GLOV® SUPERNOVA Led Patch<br>1 × Mini Controller<br>1 × Charging Case<br>1 × USB Cable<br>1 × Set of Adhesive Stickers<br>1 × User Manual | 1 × GLOV® SUPERNOVA Led Patch<br>1 × Mini Controller<br>1 × USB Cable<br>1 × Set of Adhesive Stickers<br>1 × User Manual |

### 3. Description of the Device

**Light Source:** Blue (415 nm) + Red (633 nm) LEDs

**Session Duration:** 3 minutes (auto shut-off)

**Battery Life:** 8–12 sessions per charge

**Design:** Ergonomic **star-shaped** patch for precise spot care

**Frequency:** Safe for daily use; multiple areas can be treated consecutively

### 4. Before You Start

Perform a sensitivity test before first use:

1. Place on inner forearm.
2. Run one full session (3 minutes).
3. If redness or irritation occurs within 24h, discontinue use.

### 5. How to Use SUPERNOVA

1. **Charge Fully** – Before first use, connect via USB or place in charging case and charge ~2 hours until powered.
2. **Prepare Skin** – Cleanse and dry the area you wish to treat. Ensure skin is free of oils, creams, or makeup for best contact.
3. **Apply Patch** – Attach the adhesive sticker to the back of the patch. Place directly on the desired skin area.
4. **Start Session** – Press and hold the power button for 1.5 seconds. The device will activate and run a **3-minute** automatic session.
5. **End of Session** – The device switches off automatically after 3 minutes. Remove patch and discard adhesive sticker.
6. **After Care** – Apply your preferred serum or moisturizer to complete your skincare routine.

△ **Note:** Multiple skin areas may be treated consecutively in the same routine.

### 6. Warnings

- Stop use if you experience discomfort, dizziness, or prolonged redness.
- Do not use on open wounds, ulcers, post-surgical areas, or severely irritated skin.
- Avoid use if you have sunburn, photosensitivity, inflamed skin, or are on photosensitizing medication.
- Do not look directly into the LEDs.
- Not for children. Use only under supervision if the user has reduced physical/mental ability.
- For hygiene reasons, do not share the patch.
- Keep away from direct sunlight, heat, or water.
- Use only the supplied USB cable and adapter. Ensure charging port is dry before use.
- Do not use while charging.
- Do not disassemble the device — it contains no user-serviceable parts.
- Maximum session time is 3 minutes (preset).

- Dispose of responsibly (see Disposal section).

## 7. Potential Benefits

Regular use of SUPERNOVA may help to:

- ✦ Purify the skin by reducing acne-causing bacteria
- ✦ Reduce the appearance of redness & inflammation
- ✦ Accelerate skin recovery post-breakout
- ✦ Support smoother, clearer-looking skin
- ✦ Complement daily skincare routines

### 7.1 Tips for Use

- Use on clean, dry skin.
- Apply serum or moisturizer after sessions.
- Safe for daily use; different areas can be treated one after another.

## 8. Using & Caring for SUPERNOVA

### 8.1 General Use

- The device works at safe, non-invasive light levels.
- Shuts off automatically after 3 minutes.

### 8.2 Cleaning

- Wipe with a damp cloth or 70% alcohol wipe.
- Do not use solvents or harsh cleaners.
- Never clean while charging.

### 8.3 Storing

- Store in the charging case.
- Avoid bending, crushing, or extreme temperatures.

### 8.4 Troubleshooting

- Device won't turn on: Recharge for 2 hours, then restart.
- Device won't turn off / frozen: Hold power button 10 seconds to reset.
- Stops during session: Recharge fully, clean contacts, retry.
- Won't charge: Use original cable, ensure port is dry, use 5V/2A adapter.

## 9. Specifications

### LED Technology

- Wavelengths: 415 nm (blue), 633 nm (red)
- Max. power density:  $\leq 5.5 \text{ mW/cm}^2$

- Session: 3 minutes (auto shut-off)

### **Battery**

- Type: Rechargeable Li-ion
- Capacity: 65 mAh
- Battery life: 8–12 sessions per charge
- Charging time: ~2 hours
- Power: 5V / 2A
- Standby: Up to 30 days

### **Physical Properties**

- Materials: PC + ABS, body-safe silicone
- Color: White/Pink
- Weight: ~18 g per patch
- Dimensions: ~35 × 35 × 2 mm

### **Environmental Conditions**

- Storage temp: 0°C to 40°C
- Storage humidity: 30% – 85%

## **10. Travel and International Use**

- Lightweight & compact – ideal for travel.
- Worldwide voltage (AC100–240V, 50/60Hz).
- Confirm airline rules before onboard use.

## **11. Limited Warranty – 2 Years**

Phenicoptere Sp. z o.o. warrants SUPERNOVA for **2 years** from purchase against defects in materials or workmanship under normal use.

### **Does not cover:**

- Cosmetic wear and tear
- Damage from misuse, neglect, or improper storage
- Unauthorized repair or tampering

⚠ Opening or disassembling voids the warranty.

If a defect occurs, Phenicoptere Sp. z o.o. may repair or replace the device. Proof of purchase required. Shipping costs non-refundable.

This warranty is in addition to your statutory consumer rights.

### **Claims:**

 [contact@phenicoptere.com](mailto:contact@phenicoptere.com)

 [www.glov.co](http://www.glov.co)

## 12. Disposal

The crossed-out bin symbol means this product must not be disposed of with household waste. Take to an official collection point for recycling.

### Battery Disposal

⚠ **CAUTION:** Opening the device voids warranty. Remove the battery only when disposing of the device.

- Contains a lithium-ion battery — do not discard with household waste.
- Switch off and disconnect before opening housing.
- Remove and recycle the battery at an official collection point.
- Wear gloves when handling.

## 13. Symbols

⚠ – Warnings / safety instructions

📖 – Refer to manual

♻ – Electronic recycling

🏭 – Manufacturer

LOT – Batch code

☔ – Keep dry

EC REP – EU Representative

## 14. Contact Information

**Manufacturer:** Produced for Phenicoptere Sp. z o.o.

**Address:** Chocimska 14a/101, 00-791 Warsaw, Poland

**Tel:** +48 577 111 675

**EU Representative:** Phenicoptere Sp. z o.o.

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## 3. Description de l'appareil

Source lumineuse : LED bleue (415 nm) + LED rouge (633 nm)

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